

## ELECTRO-CAPS

These caps have 20 tin electrodes in the standard 10/20 placements. They come with cap straps that attach to our body harness. The cap straps can be replaced with chin straps upon request.

<b>E1-L</b> .....	Electro-Cap - Large (58-62 cm) (Blue).....	\$313.00
<b>E1-L/M</b> .....	Electro-Cap - Large/Medium (56-60 cm) (Blue & Red).....	\$313.00
<b>E1-M</b> .....	Electro-Cap - Medium (54-58 cm) (Red) .....	\$313.00
<b>E1-M/SM</b> .....	Electro-Cap - Medium/Small (52-56 cm) (Red & Yellow).....	\$313.00
<b>E1-SM</b> .....	Electro-Cap - Small (50-54 cm) (Yellow).....	\$313.00
<b>E1-SM/XSM</b> .....	Electro-Cap - Small/Extra Small (48-52 cm) (Yellow & Green) .....	\$313.00
<b>E1-XSM</b> .....	Electro-Cap - Extra Small (46-50 cm) (Green) .....	\$313.00
<b>I1-1</b> .....	Infa-Cap I (42-46 cm) (Blue) (10/20 System).....	\$307.25
<b>I1-2</b> .....	Infa-Cap II (38-42 cm) (Pink) (Modified System) .....	\$283.00
<b>I1-3</b> .....	Infa-Cap III (34-38 cm) (Brown) (Modified System) .....	\$283.00

(Smaller size Infa-Caps are available)

## ELECTRO-CAP with EAR DROPS

These caps have the same electrode layout as the Electro-Caps, but also have dropped wires to accommodate ear electrodes (item E5-9S DROPS). They also come with cap straps, which can be replaced with chin straps upon request.

<b>E1-LEX-(size)</b> .....	Electro-Cap with ear drops (specify size).....	\$342.40
----------------------------	--	----------

## SURGICAL ELECTRO-CAP

The Surgical Cap is a strapless Electro-Cap designed to be used for intraoperative EEG monitoring. The cap fits the head securely, thereby eliminating the need for cap straps or chin straps.

<b>CS2000-(size)</b> .....	Surgical Cap (specify size).....	\$313.00
----------------------------	----------------------------------	----------

## SURGICAL ELECTRO-CAP with EAR DROPS

These caps have the same electrode layout as the surgical Electro-Caps, but also have dropped wires to accommodate ear electrodes (item E5-9S DROPS).

<b>CS2005-(size)</b> .....	Surgical Cap with ear drops (specify size) .....	\$342.40
----------------------------	--	----------

## CADWELL EASY 2 CAPS

These caps have the same electrode layout as the Electro-Caps, but have a connector that allows them to connect directly to the Cadwell Easy2 and Cadwell Easy3 equipment. They also come with cap straps, which can be replaced with chin straps upon request.

<b>E1-Easy-(size)</b> .....	Cadwell Easy2 Cap (specify size).....	\$362.90
-----------------------------	---------------------------------------	----------

## SPECIAL ORDER ELECTRO-CAPS

Special Order Electro-Caps, in standard or special sizes, with customer specified electrode locations, with electrodes made of either tin or silver/silver chloride, can be provided. Special Order Electro-Cap Forms are available. Special Order Electro-Caps are priced individually, based on the total number and placement of electrodes, the metal used for the electrodes, and any other special features that are to be incorporated into the cap design.

# PRICE LIST

May 1, 2024



### TERMS

FOB Eaton, Ohio; Net 20 days: 1 1/2% service charge per month on past due accounts.  
MINIMUM ORDER \$20.00

### LIMITED WARRANTY

Any item that fails to function properly within ninety (90) days of receipt due to defects in materials or workmanship will be repaired or replaced, at Electro-Cap International's discretion.

All prices and specifications are subject to change without notice.  
All prices are in U.S. dollars.



ISO13485:2016

## Electro-Cap International, Inc.

1011 W. Lexington Rd • P.O. Box 87 • Eaton, Ohio 45320  
800-527-2193 • 937-456-6099 • FAX 937-456-7323  
eci@electro-cap.com www.electro-cap.com

## SYSTEM PRICES

### SYS-I \$578.00

#### ELECTRO-CAP SYSTEM I

- Medium Cap (1)
- Electrode Board Adapter (1)
- Body Harness (m) (1)
- Quick Insert Electrode (1)
- Ear Electrodes (1 pr.)
- Disposable Sponge Disks (1 pk. 100)
- Needle/Syringe Kit (1)
- Electro Gel (1-16 oz.)
- Special Head Measuring Tape (1)
- Ivory® Liquid (1)
- Instruction Manual (1)
- Instructional Video (1)

### SYS-II \$944.00

#### ELECTRO-CAP SYSTEM II

- Large Cap (1)
- Medium Cap (1)
- Electrode Board Adapter (1)
- Body Harness (m) (1)
- Quick Insert Electrode (2)
- Ear Electrodes (2 pr.)
- Disposable Sponge Disks (1 pk. 100)
- Needle/Syringe Kit (2)
- Electro Gel (2-16 oz.)
- Special Head Measuring Tape (1)
- Ivory® Liquid (1)
- Instruction Manual (1)
- Instructional Video (1)

### SYS-III \$1,672.00

#### ELECTRO-CAP SYSTEM III

- Large Cap (1)
- Medium Cap (1)
- Small Cap (1)
- Extra Small Cap (1)
- Electrode Board Adapter (1)
- Body Harness (1 m; 1 sm)
- Quick Insert Electrode (4)
- Ear Electrodes (2 pr.)
- Disposable Sponge Disks (2 pk. 100)
- Needle/Syringe Kit (4)
- Electro-Gel (3-16 oz.)
- Special Head Measuring Tape (1)
- Ivory® Liquid (1)
- Instruction Manual (1)
- Instructional Video (1)

### INFA \$902.00

#### INFANT ELECTRO-CAP SYSTEM\*

- Infa-Cap I (42-46 cm ) (Blue) (10/20 System) (1)
- Infa-Cap II (38-42 cm) (Pink) (Modified System) (1)
- Infa-Cap III (34-38 cm) (Brown) (Modified System) (1)
- Infa-Body Harness (1)
- Infa-Head Measuring Tape (1)
- Skin Prep Pads (pkg/20)

Infa-Cap I has the standard international 10-20 electrode placements. Infa-Caps II and III have double interelectrode distances (electrodes are placed at F1, F2, C3, C4, O1, O2, T3, T4, CZ, and Ground).

\*Infa-Cap System requires standard E2S Electrode Board Adapter, E4S Quick Insert Electrodes, E5-6S Ear Electrodes, E6 Disposable Sponge Disks, E7 Needle/Syringe Kits, E9 Electro-Gel in addition to items listed in Infa-Cap System.

## INDIVIDUAL ITEM PRICE LIST

<b>E2S</b> .....	Electrode Board Adapter (E.B.A.) (Sockets).....	\$126.50
<b>E3-L</b> .....	Body Harness (Large) (44-52 in) .....	\$ 34.35
<b>E3-M</b> .....	Body Harness (Medium) (30-44 in).....	\$ 34.35
<b>E3-SM</b> .....	Body Harness (Small) (21-30 in).....	\$ 34.35
<b>E4S</b> .....	Quick Insert Electrode (Sockets) .....	\$ 10.65
<b>E5-9S</b> .....	Ear Electrodes, 9mm (pr.) (Sockets).....	\$ 25.30
<b>E5-6S</b> .....	Ear Electrodes, 6mm (pr.) (Sockets).....	\$ 25.30
<b>E5-9S Drops</b>	Ear Electrodes 3.5" Drops (pr.) (Sockets) .....	\$ 25.30
<b>E6</b> .....	Disposable Sponge Disks (pkg. of 100).....	\$ 24.00
<b>E7</b> .....	Needle/Syringe Kit (2-16g blunted needles; 1 syringe) (sterile).....	\$ 3.50
<b>E7B</b> .....	Syringe (5cc) (Pkg. of 125) (sterile) .....	\$120.35
** <b>E8B</b> .....	Special Blunted Needles (16g) (pkg. of 100) (sterile) .....	\$ 96.40
** <b>E8C</b> .....	Blunt Plastic Cannula (pkg. of 100) (sterile).....	\$ 96.40
<b>E9</b> .....	Electro-Gel (16 oz.).....	\$ 14.40
<b>E10</b> .....	Electro-Gel (32 oz.).....	\$ 27.70
<b>E11</b> .....	Electro-Gel (128 oz.).....	\$ 92.65
<b>E12</b> .....	Special Head Measuring Tape .....	\$ 14.90
<b>E13</b> .....	Cap Straps (Specify Large, Medium, Small, Extra Small) .....	\$ 18.10
<b>E16</b> .....	Ivory® Detergent (24 oz.).....	\$ 6.90
<b>E17S</b> .....	Bio-Potential Electrode (Sockets) .....	\$ 12.50
<b>E19</b> .....	Spare Electrode Disk .....	\$ 1.60
<b>E20</b> .....	Spare Electrode Mount .....	\$ 2.00
<b>E21-9S</b> .....	Disk Electrode, 9mm cup (Sockets).....	\$ 8.85
<b>E21-6S</b> .....	Disk Electrode, 6mm cup (Sockets).....	\$ 8.85
<b>E22</b> .....	Extra Electrodes added to cap (ea).....	\$ 17.85
<b>E23-5</b> .....	Chin Straps 5" (pr.) .....	\$ 13.85
<b>E23-7</b> .....	Chin Straps 7" (pr.) .....	\$ 13.85
<b>E24</b> .....	Chin Strap Pads (pkg. of 100).....	\$ 24.00
<b>E25</b> .....	Adhesive Electrode Pads (pkg. of 100).....	\$ 24.00
<b>E26</b> .....	Orange Sticks (pkg. of 100) .....	\$ 10.65
<b>E28S</b> .....	Alligator Clip Adapter (pkg. of 2) (Sockets).....	\$ 25.30
<b>E30B</b> .....	Syringe (10cc) (Sterile) (pkg. of 200) .....	\$169.35
<b>E40</b> .....	Converter - Pin/Socket.....	\$ 2.95
<b>I2</b> .....	Infa-Body Harness .....	\$ 19.20
<b>I3</b> .....	Infa Head Measuring Tape.....	\$ 13.85
<b>I4</b> .....	Skin Prep Pads (pkg. of 20).....	\$ 6.25
<b>INS-1</b> .....	Instruction Manual.....	\$ 12.25
<b>VID-1</b> .....	Instructional DVD .....	\$ 18.10
<b>R1</b> .....	Repair Electro-Cap (standard 10/20 placements).....	\$118.50
<b>R2</b> .....	Repair Electrode Board Adapter (25 pin connector/20 wires).....	\$ 48.20

\*\* These items are not CE marked.

#### NuPREP, TEN20 PASTE, METRICIDE

<b>R5-4J</b> .....	Ten20 Paste (4 oz. jar) (pkg./3) .....	(8 or more \$22.75).....	\$23.30
<b>R5-4T</b> .....	Ten20 Paste (4 oz. tube) (pkg./3) .....	(8 or more \$22.75).....	\$23.30
<b>R5-8J</b> .....	Ten20 Paste (8 oz. jar) (pkg./3) .....	(8 or more \$43.90).....	\$44.75
<b>R6</b> .....	Nuprep (25 gm tubes) (pkg./6).....	(4 boxes or more \$13.90/box).....	\$14.35
<b>R7</b> .....	Nuprep (4 oz. tubes) (pkg./3).....	(8 or more \$25.50).....	\$26.20
<b>R8-2800</b> .....	MetriCide 28 (30 day gal.).....		\$28.60
<b>R8-2805</b> .....	MetriCide 28 (30 day qt.) .....		\$11.80